



SURGICAL SUTURES





Suture materials

ABSORBABLE SUTURES

Loss of tensile strength over time:




SYNTHETIC ABSORBABLE SURGICAL SUTURE
LUXCRYL 910
POLYGLACTINE 910 BRAIDED




SYNTHETIC ABSORBABLE SURGICAL SUTURE
LUXCRYL PGA
POLYGLYCOLIC ACID BRAIDED




SYNTHETIC ABSORBABLE SURGICAL SUTURE
LUXCRYL PGA RAPID
POLYGLYCOLIC ACID BRAIDED




SYNTHETIC ABSORBABLE SURGICAL SUTURE
LUXCRYL PDO
POLYDIOXANONE MONOFILAMENT



SYNTHETIC ABSORBABLE SURGICAL SUTURE
LUXCRYL MONOFAST
POLYGLCAPRONE 25



ABSORBABLE SURGICAL SUTURE
CATGUT PLAIN




ABSORBABLE SURGICAL SUTURE
CATGUT CHROM



NON-ABSORBABLE SUTURES




NONABSORBABLE SURGICAL SUTURE
LUXYLENE
POLYPROPYLENE MONOFILAMENT




NONABSORBABLE SURGICAL SUTURE
LUXPET
POLYESTER BRAIDED




NONABSORBABLE SURGICAL SUTURE
LUXAMID
NYLON BLUE



NONABSORBABLE SURGICAL SUTURE
SILK BRAIDED
BLACK




NONABSORBABLE SURGICAL SUTURE
LUXAMID
NYLON BLACK



NONABSORBABLE SURGICAL SUTURE
SILK VIRGIN
BLUE



NONABSORBABLE SURGICAL SUTURE
SUPRAMID
POLYAMID BLACK




NONABSORBABLE SURGICAL SUTURE
STEEL
STAINLESS STEEL MONOFILAMENT

ORDER NUMBER SPECIFICATION – FOIL PACKS



PL 30 HR 22 - X



SUTURE MATERIAL

PL	Luxcryl 910 polyglactine 910
PG	Luxcryl PGA polyglycolic acid
PGR	Luxcryl PGA Rapid polyglycolic acid
PD	Luxcryl PDO polydioxanone
PC	Luxcryl Monofast polyglecaprone 25
CP	Catgut Plain
CC	Catgut Chrom
PP	Luxylene polypropylene
NY	Luxamid nylon
SB	Supramid polyamid
SI	Silk black
SV	Silk virgin
PE	Luxpet polyester
ST	Steel

NEEDLE LENGTH

mm

THREAD SIZE

USP

OTHER SPECIFICATION

D	double needle
L	loop
P	premium needle
U	undyed suture
UP	undyed suture, premium needle

NEEDLE TYPE

●	HR	1/2 circle • round bodied • taper point
	DR	3/8 circle • round bodied • taper point
	VR	1/4 circle • round bodied • taper point
	FR	5/8 circle • round bodied • taper point
	GR	straight • round bodied • taper point
	JR	J-needle • round bodied • taper point
⊙	HRT	1/2 circle • round bodied • tapercut point
	DRT	3/8 circle • round bodied • tapercut point
	JRT	J-needle • round bodied • tapercut point
▼	HS	1/2 circle • reverse cutting
	DS	3/8 circle • reverse cutting
▲	NH	1/2 circle • cutting
	NS	3/8 circle • cutting
	GS	straight • cutting
◄	HSP	1/2 circle • side cutting spatula
	DSP	3/8 circle • side cutting spatula
	VSP	1/4 circle • side cutting spatula
	GSP	straight • side cutting spatula
●	I CC	3/8 circle • soft cutting point • coronary
●	I DRB	3/8 circle • blunt needle
○	I WN	without needle

NEEDLE TYPES




5/8 circle




1/2 circle




3/8 circle




1/4 circle



J-needle




straight




ROUND BODIED TAPER POINT

A gradual taper to a round bodied needle.




SIDE CUTTING SPATULA NEEDLES

Needle with cutting edges on the two sides.




ROUND BODIED TAPERCUT POINT

Balanced, cutting point-tapering into a round bodied needle.




SOFT CUTTING POINT - CORONARY

A minute point followed by a round body and a rectangular main body for better penetration in cardiovascular surgery.




REVERSE CUTTING

A triangular body offers a third cutting edge on the outer convex curvature.




CONVENTIONAL CUTTING

A triangular body offers a third cutting edge on the inside concave curvature.



PREMIUM NEEDLE

Precision point with extra sharpening.



BLUNT NEEDLE

Cylindrical bodied needle with a rounded point. For soft tissue such as in liver and kidneys. Reduced risk of needle stick injuries.

SYNTHETIC ABSORBABLE SURGICAL SUTURE

LUXCRYL 910

POLYGLACTINE 910 BRAIDED

TYPE	Braided and coated multifilament	
COMPOSITION	Polyglactine 910, a copolymer made of 90% glycolide and 10%L-lactide	
COATING	Poly-glycolide-co-L-lactide and calcium stearate (<1%)	
TISSUE REACTION	Minimal	
ABSORPTION	The hydrolytic action by which the material is broken down results in total absorption between 56 and 70 days. Between 40% and 50% of tensile strength remain after 21 days.	
TENSILE STRENGTH	Post implantation	Approximate % original strength remaining
	14 days	75 %
	21 days	50 %
	28 days	20 %
SUTURE COLOUR	Violet	
RANGE	USP 8/0 - 6	
INDICATIONS	Luxcryl 910 sutures are intended for use in general soft tissue closing and/or ligation; especially in general surgery, skin closure, gastrointestinal surgery, plastic surgery, ophthalmic surgery, gynaecology, obstetrics, orthopaedics, urology	



Without needles

NEEDLE	EP	USP	LENGTH	SUTURE COLOUR	OTHER INFORMATION	ORDER NUMBER
	1.5	4/0	5 x 70 cm	violet		PL40WN70-5
	2	3/0	5 x 70 cm	violet		PL30WN70-5
	3	2/0	5 x 70 cm	violet		PL20WN70-5
	3.5	0	5 x 70 cm	violet		PL00WN70-5
	4	1	5 x 70 cm	violet		PL01WN70-5
	1.5	4/0	150 cm	violet		PL40WN150
	2	3/0	150 cm	violet		PL30WN150
	3	2/0	150 cm	violet		PL20WN150
	3.5	0	150 cm	violet		PL00WN150
	4	1	150 cm	violet		PL01WN150
	5	2	150 cm	violet		PL02WN150

Without needles – ligapack

NEEDLE	EP	USP	LENGTH	SUTURE COLOUR	OTHER INFORMATION	ORDER NUMBER
	2	3/0	250 cm	violet		PL30WN250
	3	2/0	250 cm	violet		PL20WN250
	3.5	0	250 cm	violet		PL00WN250
	4	1	250 cm	violet		PL01WN250
	5	2	250 cm	violet		PL02WN250



Round bodied needle – taper point – 1/2 circle

NEEDLE	EP	USP	LENGTH	SUTURE COLOUR	OTHER INFORMATION	ORDER NUMBER
HR13 13 mm 	0.7	6/0	75 cm	violet		PL60HR13
	1	5/0	75 cm	violet		PL50HR13
	1.5	4/0	75 cm	violet		PL40HR13
HR17 17 mm 	1	5/0	75 cm	violet		PL50HR17
	1.5	4/0	75 cm	violet		PL40HR17
	2	3/0	75 cm	violet		PL30HR17
HR20 20 mm 	1.5	4/0	75 cm	violet		PL40HR20
HR22 22 mm 	1.5	4/0	75 cm	violet		PL40HR22
	2	3/0	75 cm	violet		PL30HR22
	3	2/0	75 cm	violet		PL20HR22



Round bodied needle – taper point – 1/2 circle

NEEDLE	EP	USP	LENGTH	SUTURE COLOUR	OTHER INFORMATION	ORDER NUMBER
HR26 26 mm 	1.5	4/0	75 cm	violet		PL40HR26
	2	3/0	75 cm	violet		PL30HR26
	3	2/0	75 cm	violet		PL20HR26
	3	2/0	90 cm	violet		PL20HR26-90
	3.5	0	75 cm	violet		PL00HR26
	4	1	75 cm	violet		PL01HR26
HR30 30 mm 	2	3/0	75 cm	violet		PL30HR30
	3	2/0	75 cm	violet		PL20HR30
	3	2/0	90 cm	violet		PL20HR30-90
	3.5	0	75 cm	violet		PL00HR30
	4	1	90 cm	violet		PL01HR30
HR36 36 mm 	2	3/0	75 cm	violet		PL30HR36
	3	2/0	75 cm	violet		PL20HR36
	3	2/0	90 cm	violet		PL20HR36-90
	3.5	0	90 cm	violet		PL00HR36
	4	1	90 cm	violet		PL01HR36
	5	2	90 cm	violet		PL02HR36
HR40 40 mm 	3	2/0	75 cm	violet		PL20HR40
	3	2/0	90 cm	violet		PL20HR40-90
	3.5	0	90 cm	violet		PL00HR40
	4	1	90 cm	violet		PL01HR40
	4	1	75 cm	violet		PL01HR40-75
	5	2	90 cm	violet		PL02HR40
HR45 45 mm 	4	1	90 cm	violet		PL01HR45
	5	2	90 cm	violet		PL02HR45
HR48 48 mm 	3.5	0	90 cm	violet		PL00HR48
	4	1	90 cm	violet		PL01HR48
	5	2	90 cm	violet		PL02HR48



Round bodied needle – blunt – 3/8 circle

NEEDLE	EP	USP	LENGTH	SUTURE COLOUR	OTHER INFORMATION	ORDER NUMBER
DRB63 63 mm 	4	1	100 cm	violet		PL01DRB63



Round bodied needle – tapercut point – 1/2 circle

NEEDLE	EP	USP	LENGTH	SUTURE COLOUR	OTHER INFORMATION	ORDER NUMBER
HRT18 18 mm 	1.5	4/0	75 cm	violet		PL40HRT18



Round bodied needle – tapercut point – J-needle

NEEDLE	EP	USP	LENGTH	SUTURE COLOUR	OTHER INFORMATION	ORDER NUMBER
JRT30 30 mm 	4	1	90 cm	violet		PL01JRT30



Reverse cutting needle – 3/8 circle

NEEDLE	EP	USP	LENGTH	SUTURE COLOUR	OTHER INFORMATION	ORDER NUMBER
DS12 12 mm 	0.7	6/0	45 cm	violet		PL60DS12
	0.7	6/0	45 cm	undyed		PL60DS12-U
	1	5/0	45 cm	violet		PL50DS12
	1	5/0	45 cm	undyed		PL50DS12-U
	1.5	4/0	45 cm	violet		PL40DS12
DS16 16 mm 	0.7	6/0	45 cm	violet		PL60DS16
	0.7	6/0	45 cm	undyed		PL60DS16-U
	1	5/0	75 cm	violet		PL50DS16
	1.5	4/0	75 cm	violet		PL40DS16
	1.5	4/0	75 cm	undyed		PL40DS16-U
	2	3/0	75 cm	violet		PL30DS16
DS19 19 mm 	1	5/0	75 cm	violet		PL50DS19
	1	5/0	75 cm	undyed		PL50DS19-U
	1.5	4/0	75 cm	violet		PL40DS19
	1.5	4/0	75 cm	undyed		PL40DS19-U
	2	3/0	75 cm	violet		PL30DS19
	2	3/0	75 cm	undyed		PL30DS19-U
DS24 24 mm 	1.5	4/0	75 cm	violet		PL40DS24
	2	3/0	75 cm	violet		PL30DS24
	2	3/0	75 cm	undyed		PL30DS24-U
	2	3/0	90 cm	violet		PL30DS24-90
	3	2/0	75 cm	violet		PL20DS24
DS26 26 mm 	3	2/0	90 cm	violet		PL20DS26-90
	3	2/0	75 cm	undyed		PL20DS26-U

Reverse cutting needle – 3/8 circle						
NEEDLE	EP	USP	LENGTH	SUTURE COLOUR	OTHER INFORMATION	ORDER NUMBER
DS30 30 mm 	2	3/0	75 cm	violet		PL30DS30
	3	2/0	75 cm	violet		PL20DS30
	3	2/0	75 cm	undyed		PL20DS30-U
	3	2/0	90 cm	violet		PL20DS30-90
	3.5	0	75 cm	violet		PL00DS30
	3.5	0	75 cm	undyed		PL00DS30-U
DS36 36 mm 	4	1	75 cm	violet		PL01DS30
	3	2/0	75 cm	violet		PL20DS36
	3	2/0	75 cm	undyed		PL20DS36-U
	3	2/0	90 cm	violet		PL20DS36-90
	3.5	0	75 cm	violet		PL00DS36
DS40 40 mm 	4	1	75 cm	violet		PL01DS36
	3.5	0	75 cm	violet		PL00DS40
	3.5	0	90 cm	violet		PL00DS40-90
	4	1	75 cm	violet		PL01DS40
	4	1	90 cm	violet		PL01DS40-90

Reverse cutting needle – 1/2 circle						
NEEDLE	EP	USP	LENGTH	SUTURE COLOUR	OTHER INFORMATION	ORDER NUMBER
HS26 26 mm 	2	3/0	75 cm	violet		PL30HS26
	3	2/0	75 cm	violet		PL20HS26
HS36 36 mm 	3	2/0	75 cm	violet		PL20HS36
	3.5	0	90 cm	violet		PL00HS36
	4	1	90 cm	violet		PL01HS36
HS40 40 mm 	3.5	0	90 cm	violet		PL00HS40
	4	1	90 cm	violet		PL01HS40
HS48 48 mm 	3.5	0	90 cm	violet		PL00HS48
	4	1	90 cm	violet		PL01HS48
	5	2	90 cm	violet		PL02HS48

Reverse cutting needle – straight						
NEEDLE	EP	USP	LENGTH	SUTURE COLOUR	OTHER INFORMATION	ORDER NUMBER
GS60 60 mm 	3	2/0	75 cm	violet		PL20GS60

Sutures in cassettes



EP	USP	LENGTH	SUTURE COLOUR	OTHER INFORMATION	ORDER NUMBER
1.5	4/0	15 m	violet		PL4015
2	3/0	15 m	violet		PL3015
2	3/0	25 m	violet		PL3025
3	2/0	15 m	violet		PL3015
3	2/0	25 m	violet		PL3025
3	2/0	50 m	violet		PL3050
3.5	0	15 m	violet		PL0015
3.5	0	25 m	violet		PL0025
3.5	0	50 m	violet		PL0050
4	1	15 m	violet		PL0115
4	1	25 m	violet		PL0125
4	1	50 m	violet		PL0150
5	2	15 m	violet		PL0215
5	2	25 m	violet		PL0225
5	2	50 m	violet		PL0250
5	2	100 m	violet		PL02100
6	3&4	15 m	violet		PL0315
6	3&4	25 m	violet		PL0325
6	3&4	50 m	violet		PL0350
7	5	15 m	violet		PL0515
7	5	25 m	violet		PL0525
7	5	50 m	violet		PL0550
8	6	15 m	violet		PL0615
8	6	25 m	violet		PL0625
8	6	50 m	violet		PL0650

SYNTHETIC ABSORBABLE SURGICAL SUTURE

LUXCRYL PGA

POLYGLYCOLIC ACID BRAIDED

TYPE	Braided and coated multifilament	
COMPOSITION	Polyglycolic acid	
COATING	Polycaprolactone and calcium stearate (1%)	
TISSUE REACTION	Minimal	
ABSORPTION	The hydrolytic action by which the material is broken down results in total absorption in approximately 60 to 90 days. Approximately 50% of tensile strength remain after 21 days.	
TENSILE STRENGTH	Post implantation	Approximate % original strength remaining
	14 days	70 %
	21 days	50 %
	28 days	20 %
SUTURE COLOUR	Violet, beige - undyed	
RANGE	USP 10/0 - 6	
INDICATIONS	Luxcryl PGA sutures are intended for use in general soft tissue closing and/or ligation; especially in general surgery, skin closure, gastrointestinal surgery, gynaecology, obstetrics, plastic surgery, urology, ophthalmic surgery, orthopaedics	



Without needles – ligapack

NEEDLE	EP	USP	LENGTH	SUTURE COLOUR	OTHER INFORMATION	ORDER NUMBER
	2	3/0	250 cm	violet		PG30WN250
	3	2/0	250 cm	violet		PG20WN250
	3.5	0	250 cm	violet		PG00WN250
	4	1	250 cm	violet		PG01WN250



Round bodied needle – taper point – 1/2 circle

NEEDLE	EP	USP	LENGTH	SUTURE COLOUR	OTHER INFORMATION	ORDER NUMBER
HR13 13 mm 	0.7	6/0	75 cm	violet		PG60HR13
HR17 17 mm 	1	5/0	75 cm	violet		PG50HR17
	1.5	4/0	75 cm	violet		PG40HR17
HR22 22 mm 	1.5	4/0	75 cm	violet		PG40HR22
	2	3/0	75 cm	violet		PG30HR22
	3	2/0	75 cm	violet		PG20HR22
HR26 26 mm 	1.5	4/0	75 cm	violet		PG40HR26
	2	3/0	75 cm	violet		PG30HR26
	3	2/0	75 cm	violet		PG20HR26
HR30 30 mm 	2	3/0	75 cm	violet		PG30HR30
	3	2/0	75 cm	violet		PG20HR30
	3.5	0	75 cm	violet		PG00HR30
HR36 36 mm 	2	3/0	75 cm	violet		PG30HR36
	3	2/0	75 cm	violet		PG20HR36
	3.5	0	90 cm	violet		PG00HR36
	4	1	90 cm	violet		PG01HR36
HR40 40 mm 	3.5	0	90 cm	violet		PG00HR40
	4	1	90 cm	violet		PG01HR40
	5	2	90 cm	violet		PG02HR40
HR48 48 mm 	3.5	0	90 cm	violet		PG00HR48
	4	1	90 cm	violet		PG01HR48
	5	2	90 cm	violet		PG02HR48

Reverse cutting needle – 3/8 circle							
NEEDLE		EP	USP	LENGTH	SUTURE COLOUR	OTHER INFORMATION	ORDER NUMBER
DS12	12 mm		0.7	6/0	45 cm	violet	PG60DS12
			1	5/0	45 cm	violet	PG50DS12
DS16	16 mm		0.7	6/0	75 cm	violet	PG60DS16
			1	5/0	75 cm	violet	PG50DS16
			1.5	4/0	75 cm	violet	PG40DS16
			2	3/0	75 cm	violet	PG30DS16
DS19	19 mm		1	5/0	75 cm	violet	PG50DS19
			1.5	4/0	75 cm	violet	PG40DS19
			2	3/0	75 cm	violet	PG30DS19
			3	2/0	75 cm	violet	PG20DS19
DS24	24 mm		2	3/0	75 cm	violet	PG30DS24
			3	2/0	75 cm	violet	PG20DS24
DS30	30 mm		2	3/0	75 cm	violet	PG30DS30
			3	2/0	75 cm	violet	PG20DS30
			3.5	0	75 cm	violet	PG00DS30
			4	1	75 cm	violet	PG01DS30
DS40	40 mm		4	1	75 cm	violet	PG01DS40

Reverse cutting needle – 1/2 circle							
NEEDLE		EP	USP	LENGTH	SUTURE COLOUR	OTHER INFORMATION	ORDER NUMBER
HS22	22 mm		1.5	4/0	75 cm	violet	PG40HS22
HS26	26 mm		2	3/0	75 cm	violet	PG30HS26
			3	2/0	75 cm	violet	PG20HS26
HS36	36 mm		3.5	0	90 cm	violet	PG00HS36
HS40	40 mm		3.5	0	90 cm	violet	PG00HS40
			4	1	90 cm	violet	PG01HS40
			5	2	90 cm	violet	PG02HS40
HS48	48 mm		5	2	90 cm	violet	PG02HS48

Sutures in cassettes



EP	USP	LENGTH	SUTURE COLOUR	OTHER INFORMATION	ORDER NUMBER
1.5	4/0	15 m	violet		PG4015
2	3/0	15 m	violet		PG3015
2	3/0	25 m	violet		PG3025
3	2/0	15 m	violet		PG2015
3	2/0	25 m	violet		PG2025
3	2/0	50 m	violet		PG2050
3.5	0	15 m	violet		PG0015
3.5	0	25 m	violet		PG0025
3.5	0	50 m	violet		PG0050
4	1	15 m	violet		PG0115
4	1	25 m	violet		PG0125
4	1	50 m	violet		PG0150
5	2	15 m	violet		PG0215
5	2	25 m	violet		PG0225
5	2	50 m	violet		PG0250
5	2	100 m	violet		PG02100
6	3&4	15 m	violet		PG0315
6	3&4	25 m	violet		PG0325
6	3&4	50 m	violet		PG0350
6	3&4	100 m	violet		PG03100
7	5	15 m	violet		PG0515
7	5	25 m	violet		PG0525
7	5	50 m	violet		PG0550
8	6	15 m	violet		PG0615
8	6	25 m	violet		PG0625
8	6	50 m	violet		PG0650

Sutures in cassettes – with storage solution



EP	USP	LENGTH	SUTURE COLOUR	OTHER INFORMATION	ORDER NUMBER
1.5	4/0	15 m	violet	with storage solution	PGL4015
1.5	4/0	25 m	violet	with storage solution	PGL4025
2	3/0	15 m	violet	with storage solution	PGL3015
2	3/0	25 m	violet	with storage solution	PGL3025
3	2/0	15 m	violet	with storage solution	PGL2015
3	2/0	25 m	violet	with storage solution	PGL2025
3	2/0	50 m	violet	with storage solution	PGL2050
3.5	0	15 m	violet	with storage solution	PGL0015
3.5	0	25 m	violet	with storage solution	PGL0025
3.5	0	50 m	violet	with storage solution	PGL0050
4	1	15 m	violet	with storage solution	PGL0115
4	1	25 m	violet	with storage solution	PGL0125
4	1	50 m	violet	with storage solution	PGL0150
5	2	15 m	violet	with storage solution	PGL0215
5	2	25 m	violet	with storage solution	PGL0225
5	2	50 m	violet	with storage solution	PGL0250
5	2	100 m	violet	with storage solution	PGL02100
6	3&4	15 m	violet	with storage solution	PGL0315
6	3&4	25 m	violet	with storage solution	PGL0325
6	3&4	50 m	violet	with storage solution	PGL0350
6	3&4	100 m	violet	with storage solution	PGL03100
7	5	15 m	violet	with storage solution	PGL0515
7	5	25 m	violet	with storage solution	PGL0525
7	5	50 m	violet	with storage solution	PGL0550
8	6	15 m	violet	with storage solution	PGL0615
8	6	25 m	violet	with storage solution	PGL0625
8	6	50 m	violet	with storage solution	PGL0650

SYNTHETIC ABSORBABLE SURGICAL SUTURE

LUXCRYL PGA RAPID

POLYGLYCOLIC ACID BRAIDED

TYPE	Braided and coated multifilament	
COMPOSITION	Polyglycolic acid	
COATING	Polycaprolactone and calcium stearate (1%)	
TISSUE REACTION	Minimal	
ABSORPTION	The hydrolytic action by which the material is broken down results in total absorption in approximately 42 days. Approximately 50% of tensile strength remain after 7 days.	
TENSILE STRENGTH	Post implantation 7 days	Approximate % original strength remaining 50 %
SUTURE COLOUR	Beige - undyed, violet	
RANGE	USP 6/0 - 2	
INDICATIONS	Luxcryl PGA Rapid sutures are intended for use in superficial soft tissue closing of skin and mucosa only; especially in plastic surgery, obstetrics (episiotomy repair), gynaecology, urology (circumcision), ophthalmic surgery (conjunctiva), pediatric surgery	





Round bodied needle – taper point – 1/2 circle

NEEDLE			EP	USP	LENGTH	SUTURE COLOUR	OTHER INFORMATION	ORDER NUMBER
HR13	13 mm		0.7	6/0	75 cm	undyed		PGR60HR13
HR17	17 mm		1	5/0	75 cm	undyed		PGR50HR17
			1.5	4/0	75 cm	undyed		PGR40HR17
HR22	22 mm		1.5	4/0	75 cm	undyed		PGR40HR22
			2	3/0	75 cm	undyed		PGR30HR22
HR26	26 mm		2	3/0	75 cm	undyed		PGR30HR26
			3	2/0	75 cm	undyed		PGR20HR26
HR30	30 mm		3	2/0	75 cm	undyed		PGR20HR30
			3.5	0	75 cm	undyed		PGR00HR30
HR36	36 mm		2	3/0	90 cm	undyed		PGR30HR36-90
			3	2/0	90 cm	undyed		PGR20HR36-90
			3.5	0	75 cm	undyed		PGR00HR36
			4	1	75 cm	undyed		PGR01HR36
HR40	40 mm		3.5	0	75 cm	undyed		PGR00HR40
			4	1	75 cm	undyed		PGR01HR40
			5	2	75 cm	undyed		PGR02HR40
HR48	48 mm		4	1	75 cm	undyed		PGR01HR48
			5	2	75 cm	undyed		PGR02HR48



Reverse cutting needle – 3/8 circle

NEEDLE			EP	USP	LENGTH	SUTURE COLOUR	OTHER INFORMATION	ORDER NUMBER
DS12	12 mm		0.7	6/0	75 cm	undyed		PGR60DS12
			1	5/0	75 cm	undyed		PGR50DS12
DS16	16 mm		1	5/0	75 cm	undyed		PGR50DS16
			1.5	4/0	75 cm	undyed		PGR40DS16
			2	3/0	75 cm	undyed		PGR30DS16
DS19	19 mm		1	5/0	75 cm	undyed		PGR50DS19
			1.5	4/0	75 cm	undyed		PGR40DS19
			2	3/0	75 cm	undyed		PGR30DS19
DS24	24 mm		2	3/0	75 cm	undyed		PGR30DS24
			3	2/0	75 cm	undyed		PGR20DS24
DS30	30 mm		3	2/0	90 cm	undyed		PGR20DS30-90



Round bodied needle – tapercut point – 1/2 circle

NEEDLE	EP	USP	LENGTH	SUTURE COLOUR	OTHER INFORMATION	ORDER NUMBER
HRT36 36 mm 	3	2/0	90 cm	undyed		PGR20HRT36
	3.5	0	90 cm	undyed		PGR00HRT36
	4	1	90 cm	undyed		PGR01HRT36



Reverse cutting needle – 1/2 circle

NEEDLE	EP	USP	LENGTH	SUTURE COLOUR	OTHER INFORMATION	ORDER NUMBER
HS40 40 mm 	4	1	75 cm	undyed		PGR01HS40



Reverse cutting needle – straight

NEEDLE	EP	USP	LENGTH	SUTURE COLOUR	OTHER INFORMATION	ORDER NUMBER
GS60 60 mm 	2	3/0	75 cm	undyed		PGR30GS60

LUXCRYL PDO

SYNTHETIC ABSORBABLE SURGICAL SUTURE

POLYDIOXANONE MONOFILAMENT

TYPE	Monofilament	
COMPOSITION	Polydioxanone, a polymer made from polyester poly (p-dioxanone)	
COATING	None	
TISSUE REACTION	Minimal	
ABSORPTION	The hydrolytic action by which the material is broken down results in total absorption in approximately 180 to 210 days. Approximately 50% of tensile strength remain after 42 days.	
TENSILE STRENGTH	Post implantation	Approximate % original strength remaining
	28 days	70 %
	42 days	50 %
	56 days	20 %
SUTURE COLOUR	Violet	
RANGE	USP 6/0 - 3&4	
INDICATIONS	Luxcryl PDO sutures are intended for use in general soft tissue closing; especially in general surgery, gastrointestinal surgery, ophthalmic surgery, orthopaedics, abdominal wall closure, gynaecology, plastic surgery, urology	



 Round bodied needle – taper point – 1/2 circle							
NEEDLE			EP	USP	LENGTH	SUTURE COLOUR	ORDER NUMBER
HR13	13 mm		0.7	6/0	75 cm	violet	PD60HR13
HR17	17 mm		1	5/0	75 cm	violet	PD50HR17
			1.5	4/0	75 cm	violet	PD40HR17
HR22	22 mm		1.5	4/0	75 cm	violet	PD40HR22
			2	3/0	75 cm	violet	PD30HR22
HR26	26 mm		1.5	4/0	75 cm	violet	PD40HR26
			2	3/0	75 cm	violet	PD30HR26
			3	2/0	75 cm	violet	PD20HR26
HR30	30 mm		3	2/0	75 cm	violet	PD20HR30
			3.5	0	75 cm	violet	PD00HR30
HR36	36 mm		3	2/0	75 cm	violet	PD20HR36
			3.5	0	90 cm	violet	PD00HR36
			4	1	90 cm	violet	PD01HR36
HR40	40 mm		3.5	0	90 cm	violet	PD00HR40
			4	1	90 cm	violet	PD01HR40
			5	2	90 cm	violet	PD02HR40

 Round bodied needle – taper point – 1/2 circle – loop							
NEEDLE			EP	USP	LENGTH	SUTURE COLOUR	ORDER NUMBER
HR50	50 mm		4	1	150 cm	violet	PD01HR50-L

Reverse cutting needle – 3/8 circle							
NEEDLE			EP	USP	LENGTH	SUTURE COLOUR	OTHER INFORMATION
DS12	12 mm		0.7	6/0	45 cm	violet	PD60DS12
			1	5/0	45 cm	violet	PD50DS12
DS16	16 mm		0.7	6/0	45 cm	violet	PD60DS16
			1	5/0	45 cm	violet	PD50DS16
			1.5	4/0	75 cm	violet	PD40DS16
DS19	19 mm		1	5/0	45 cm	violet	PD50DS19
			1	5/0	75 cm	violet	PD50DS19-75
			1.5	4/0	75 cm	violet	PD40DS19
			2	3/0	75 cm	violet	PD30DS19
DS24	24 mm		2	3/0	75 cm	violet	PD30DS24
			3	2/0	75 cm	violet	PD20DS24
			3	2/0	90 cm	violet	PD20DS24-90
			3.5	0	75 cm	violet	PD00DS24
DS26	26 mm		2	3/0	75 cm	violet	PD30DS26
			2	3/0	90 cm	violet	PD30DS26-90
DS30	30 mm		3	2/0	75 cm	violet	PD20DS30
			3.5	0	75 cm	violet	PD00DS30
			3.5	0	90 cm	violet	PD00DS30-90
DS40	40 mm		4	1	90 cm	violet	PD01DS40

Reverse cutting needle – 1/2 circle							
NEEDLE			EP	USP	LENGTH	SUTURE COLOUR	OTHER INFORMATION
HS36	36 mm		3.5	0	75 cm	violet	PD00HS36

Sutures in cassettes



EP	USP	LENGTH	SUTURE COLOUR	OTHER INFORMATION	ORDER NUMBER
1.5	4/0	15 m	violet		PD4015
2	3/0	15 m	violet		PD3015
2	3/0	25 m	violet		PD3025
3	2/0	15 m	violet		PD2015
3	2/0	25 m	violet		PD2025
3.5	0	15 m	violet		PD0015
3.5	0	25 m	violet		PD0025
4	1	15 m	violet		PD0115
4	1	25 m	violet		PD0125
4	1	50 m	violet		PD0150
5	2	15 m	violet		PD0215
5	2	25 m	violet		PD0225
5	2	50 m	violet		PD0250
6	3&4	15 m	violet		PD0315
6	3&4	25 m	violet		PD0325
6	3&4	50 m	violet		PD0350

SYNTHETIC ABSORBABLE SURGICAL SUTURE

LUXCRYL MONOFAST

POLYGLYCAPRONE 25

TYPE	Monofilament	
COMPOSITION	Polyglycaprone 25, a copolymer made of 75 % glycolide and 25 % ε-caprolactone	
COATING	None	
TISSUE REACTION	Minimal	
ABSORPTION	The hydrolytic action by which the material is broken down results in total absorption in approximately 90 to 120 days. Approximately 60% of tensile strength remain after 7 days.	
TENSILE STRENGTH	Post implantation 7 days 14 days	Approximate % original strength remaining 60 % 30 %
SUTURE COLOUR	Violet, beige - undyed	
RANGE	USP 6/0 - 2	
INDICATIONS	Luxcryl Monofast sutures are intended for use in superficial soft tissue closing and/or ligation; especially in general surgery, skin closure, plastic surgery, gastrointestinal surgery, gynaecology, urology	





Round bodied needle – taper point – 1/2 circle

NEEDLE			EP	USP	LENGTH	SUTURE COLOUR	OTHER INFORMATION	ORDER NUMBER
HR13	13 mm		0.7	6/0	75 cm	violet		PC60HR13
HR17	17 mm		1	5/0	75 cm	violet		PC50HR17
HR22	22 mm		1.5	4/0	75 cm	violet		PC40HR22
			2	3/0	75 cm	violet		PC30HR22
HR26	26 mm		2	3/0	75 cm	violet		PC30HR26
			3	2/0	75 cm	violet		PC20HR26
HR30	30 mm		3	2/0	75 cm	violet		PC20HR30



Reverse cutting needle – 3/8 circle

NEEDLE			EP	USP	LENGTH	SUTURE COLOUR	OTHER INFORMATION	ORDER NUMBER
DS12	12 mm		0.7	6/0	45 cm	undyed		PC60DS12
			0.7	6/0	45 cm	violet		PC60DS12-V
			0.7	6/0	75 cm	violet		PC60DS12-75V
			1	5/0	45 cm	violet		PC50DS12-V
DS16	16 mm		0.7	6/0	90 cm	violet		PC60DS16-90V
			1	5/0	45 cm	undyed		PC50DS16
			1	5/0	45 cm	violet		PC50DS16-V
			1	5/0	75 cm	violet		PC50DS16-75V
			1.5	4/0	75 cm	violet		PC40DS16-V
DS19	19 mm		1	5/0	45 cm	undyed		PC50DS19
			1	5/0	75 cm	undyed		PC50DS19-75
			1	5/0	45 cm	violet		PC50DS19-V
			1	5/0	75 cm	violet		PC50DS19-75V
			1.5	4/0	75 cm	undyed		PC40DS19
			1.5	4/0	75 cm	violet		PC40DS19-V
			2	3/0	75 cm	undyed		PC30DS19
			2	3/0	75 cm	violet		PC30DS19-V
DS19	19 mm black		1	5/0	75 cm	violet		PC50DS19-75VB
DS24	24 mm		2	3/0	75 cm	undyed		PC30DS24
			2	3/0	75 cm	violet		PC30DS24-V
			3	2/0	75 cm	undyed		PC20DS24
			3	2/0	75 cm	violet		PC20DS24-V



Reverse cutting needle – 3/8 circle

NEEDLE			EP	USP	LENGTH	SUTURE COLOUR	OTHER INFORMATION	ORDER NUMBER
DS26	26 mm		2	3/0	75 cm	undyed		PC30DS26
DS30	30 mm		3.5	0	75 cm	violet		PC00DS30-V



Reverse cutting needle – straight

NEEDLE			EP	USP	LENGTH	SUTURE COLOUR	OTHER INFORMATION	ORDER NUMBER
GS60	60 mm		2	3/0	75 cm	undyed		PC30GS60

Sutures in cassettes



EP	USP	LENGTH	SUTURE COLOUR	OTHER INFORMATION	ORDER NUMBER
4	1	25 m	violet		PC0125
4	1	50 m	violet		PC0150
5	2	25 m	violet		PC0225
5	2	50 m	violet		PC0250

ABSORBABLE SURGICAL SUTURE CATGUT PLAIN

TYPE	Twisted multifilament with a monofilament appearance	
COMPOSITION	Strands of purified collagen taken from the serosal layer of selected bovines. Origin from BSE-free classified countries.	
COATING	None	
TISSUE REACTION	Moderate	
ABSORPTION	By phagocytosis. Complete mass absorption in approximately 63 days. Approximately 50% of tensile strength remain after 7 days. When used in infected tissues or in tissues with increased levels of proteolytic enzymes (stomach, cervix, vagina) it is more quickly absorbed.	
TENSILE STRENGTH	Post implantation 7 days 14 days	Approximate % original strength remaining 50 % 20 %
SUTURE COLOUR	Ivory	
RANGE	USP 5/0 - 6	
INDICATIONS	Catgut plain sutures are intended for use in general soft tissue closing and/or ligation; especially in general surgery, gastrointestinal surgery, gynaecology, obstetrics, urology, ophthalmic surgery	





Round bodied needle – taper point – 1/2 circle

NEEDLE	EP	USP	LENGTH	SUTURE COLOUR	OTHER INFORMATION	ORDER NUMBER
HR37 37 mm 	3.5	2/0	75 cm	ivory		CP20HR37



Reverse cutting needle – 3/8 circle

NEEDLE	EP	USP	LENGTH	SUTURE COLOUR	OTHER INFORMATION	ORDER NUMBER
DS16 16 mm 	1.5	5/0	75 cm	ivory		CP50DS16
	2	4/0	75 cm	ivory		CP40DS16
DS19 19 mm 	3	3/0	75 cm	ivory		CP30DS19
DS24 24 mm 	3.5	2/0	75 cm	ivory		CP20DS24

Sutures in cassettes



EP	USP	LENGTH	SUTURE COLOUR	OTHER INFORMATION	ORDER NUMBER
2	4/0	100 m	ivory		CP40100
3	3/0	100 m	ivory		CP30100
3.5	2/0	100 m	ivory		CP20100
4	0	100 m	ivory		CP00100
5	1	75 m	ivory		CP0175
6	2	50 m	ivory		CP0250
7	3	50 m	ivory		CP0350
8	4	30 m	ivory		CP0430
9	5	30 m	ivory		CP0530
10	6	30 m	ivory		CP0630

ABSORBABLE SURGICAL SUTURE

Without needles – ligapack

NEEDLE	EP	USP	LENGTH	SUTURE COLOUR	OTHER INFORMATION	ORDER NUMBER
	7	3	150 cm	brown		CC03WN150
	8	4	150 cm	brown		CC04WN150




Round bodied needle – taper point – 1/2 circle

NEEDLE	EP	USP	LENGTH	SUTURE COLOUR	OTHER INFORMATION	ORDER NUMBER
HR20 20 mm 	1.5	5/0	75 cm	brown		CC50HR20
	2	4/0	75 cm	brown		CC40HR20
	3	3/0	75 cm	brown		CC30HR20
HR26 26 mm 	2	4/0	75 cm	brown		CC40HR26
	3	3/0	75 cm	brown		CC30HR26
	3.5	2/0	75 cm	brown		CC20HR26
HR30 30 mm 	3.5	2/0	75 cm	brown		CC20HR30
HR40 40 mm 	3.5	2/0	75 cm	brown		CC20HR40
	3.5	2/0	90 cm	brown		CC20HR40-90
	4	0	75 cm	brown		CC00HR40
	4	0	90 cm	brown		CC00HR40-90
	5	1	75 cm	brown		CC01HR40
	5	1	90 cm	brown		CC01HR40-90
	6	2	75 cm	brown		CC02HR40
HR48 48 mm 	6	2	75 cm	brown		CC02HR48



Round bodied needle – taper point – 3/8 circle

NEEDLE	EP	USP	LENGTH	SUTURE COLOUR	OTHER INFORMATION	ORDER NUMBER
DR30 30 mm 	4	0	75 cm	brown		CC00DR30-75
	4	0	90 cm	brown		CC00DR30
DR45 45 mm 	5	1	90 cm	brown		CC01DR45



NEEDLE			EP	USP	LENGTH	SUTURE COLOUR	OTHER INFORMATION	ORDER NUMBER
DS12	12 mm		1	6/0	75 cm	brown		CC60DS12-75
			1.5	5/0	75 cm	brown		CC50DS12-75
DS16	16 mm		2	4/0	75 cm	brown		CC40DS16
DS19	19 mm		1.5	5/0	75 cm	brown		CC50DS19
			2	4/0	75 cm	brown		CC40DS19
			3	3/0	75 cm	brown		CC30DS19
DS26	26 mm		3	3/0	75 cm	brown		CC30DS26
			3.5	2/0	75 cm	brown		CC20DS26
DS30	30 mm		3.5	2/0	75 cm	brown		CC20DS30
DS37	37 mm		4	0	75 cm	brown		CC00DS37

Sutures in cassettes



EP	USP	LENGTH	SUTURE COLOUR	OTHER INFORMATION	ORDER NUMBER
2	4/0	100 m	brown		CC40100
2	4/0	50 m	brown		CC4050
3	3/0	100 m	brown		CC30100
3	3/0	50 m	brown		CC3050
3.5	2/0	100 m	brown		CC20100
3.5	2/0	50 m	brown		CC2050
4	0	100 m	brown		CC00100
4	0	75 m	brown		CC0075
4	0	50 m	brown		CC0050
5	1	75 m	brown		CC0175
5	1	50 m	brown		CC0150
6	2	50 m	brown		CC0250
7	3	50 m	brown		CC0350
8	4	30 m	brown		CC0430
8	4	25 m	brown		CC0430
9	5	30 m	brown		CC0530
10	6	30 m	brown		CC0630

LUXYLENE

NONABSORBABLE SURGICAL SUTURE
POLYPROPYLENE MONOFILAMENT

TYPE	Monofilament
COMPOSITION	Polypropylene, a polymer of propylene
COATING	None
TISSUE REACTION	Minimal
ABSORPTION	Non absorbable
SUTURE COLOUR	Blue
RANGE	USP 10/0 - 2
INDICATIONS	Luxylene sutures are intended for use in general soft tissue closing and/or ligation; especially in ophthalmic surgery, plastic surgery, skin closure, orthopaedics, hernias (also as a mesh), gastrointestinal surgery, general surgery, gynaecology, obstetrics



 Round bodied needle – taper point – 1/2 circle							
NEEDLE			EP	USP	LENGTH	SUTURE COLOUR	ORDER NUMBER
HR10	10 mm		0.5	7/0	45 cm	blue	PP70HR10
HR13	13 mm		0.7	6/0	75 cm	blue	PP60HR13
HR15	15 mm		1	5/0	75 cm	blue	PP50HR15
HR16	16 mm		1.5	4/0	75 cm	blue	PP40HR16
HR20	20 mm		1.5	4/0	75 cm	blue	PP40HR20
			2	3/0	75 cm	blue	PP30HR20
			1.5	4/0	75 cm	blue	PP40HR25
			2	3/0	75 cm	blue	PP30HR25
			3	2/0	75 cm	blue	PP20HR25
HR25	25 mm		3.5	0	75 cm	blue	PP00HR25
			4	1	75 cm	blue	PP01HR25
			2	3/0	75 cm	blue	PP30HR30
			3	2/0	75 cm	blue	PP20HR30
HR30	30 mm		3.5	0	75 cm	blue	PP00HR30
			3	2/0	75 cm	blue	PP20HR35
			3.5	0	75 cm	blue	PP00HR35
HR35	35 mm		4	1	75 cm	blue	PP01HR35
			3.5	0	75 cm	blue	PP00HR40
			4	1	90 cm	blue	PP01HR40
HR40	40 mm		5	2	75 cm	blue	PP02HR40

 Round bodied needle – taper point – 1/2 circle – double armed needle							
NEEDLE			EP	USP	LENGTH	SUTURE COLOUR	ORDER NUMBER
HR08	8 mm		0.4	8/0	60 cm	blue	PP80HR08-D
			0.5	7/0	60 cm	blue	PP70HR08-D
HR10	10 mm		0.5	7/0	60 cm	blue	PP70HR10-D
			0.7	6/0	60 cm	blue	PP60HR10-D
HR13	13 mm		0.7	6/0	60 cm	blue	PP60HR13-D
			1	5/0	75 cm	blue	PP50HR13-D
HR16	16 mm		1	5/0	90 cm	blue	PP50HR16-D
			1.5	4/0	90 cm	blue	PP40HR16-D
HR20	20 mm		1.5	4/0	90 cm	blue	PP40HR20-D
HR25	25 mm		1.5	4/0	90 cm	blue	PP40HR25-D
			2	3/0	90 cm	blue	PP30HR25-D
			3	2/0	90 cm	blue	PP20HR25-D



Round bodied needle – taper point – 3/8 circle

NEEDLE	EP	USP	LENGTH	SUTURE COLOUR	OTHER INFORMATION	ORDER NUMBER
DR13 13 mm 	0.7	6/0	75 cm	blue		PP60DR13



Round bodied needle – taper point – 3/8 circle – double armed needle

NEEDLE	EP	USP	LENGTH	SUTURE COLOUR	OTHER INFORMATION	ORDER NUMBER
DR08 8 mm 	0.4	8/0	45 cm	blue		PP80DR08-D
	0.5	7/0	60 cm	blue		PP70DR08-D
DR10 10 mm 	0.5	7/0	60 cm	blue		PP70DR10-D
	0.7	6/0	60 cm	blue		PP60DR10-D
DR13 13 mm 	0.7	6/0	60 cm	blue		PP60DR13-D
	1	5/0	75 cm	blue		PP50DR13-D



Reverse cutting needle – 3/8 circle

NEEDLE	EP	USP	LENGTH	SUTURE COLOUR	OTHER INFORMATION	ORDER NUMBER
DS12 12 mm 	0.7	6/0	45 cm	blue		PP60DS12
	0.7	6/0	75 cm	blue		PP60DS12-75
	1	5/0	45 cm	blue		PP50DS12
DS16 16 mm 	0.7	6/0	45 cm	blue		PP60DS16
	1	5/0	45 cm	blue		PP50DS16
	1	5/0	75 cm	blue		PP50DS16-75
	1.5	4/0	45 cm	blue		PP40DS16-45
	1.5	4/0	75 cm	blue		PP40DS16
	2	3/0	75 cm	blue		PP30DS16
DS19 19 mm 	1	5/0	45 cm	blue		PP50DS19-45
	1	5/0	75 cm	blue		PP50DS19
	1.5	4/0	45 cm	blue		PP40DS19-45
	1.5	4/0	75 cm	blue		PP40DS19
	2	3/0	75 cm	blue		PP30DS19
	3	2/0	75 cm	blue		PP20DS19
DS24 24 mm 	1.5	4/0	75 cm	blue		PP40DS24
	2	3/0	75 cm	blue		PP30DS24
	3	2/0	75 cm	blue		PP20DS24
DS30 30 mm 	2	3/0	75 cm	blue		PP30DS30
	3	2/0	75 cm	blue		PP20DS30
	3.5	0	75 cm	blue		PP00DS30

Reverse cutting needle – 3/8 circle						
NEEDLE	EP	USP	LENGTH	SUTURE COLOUR	OTHER INFORMATION	ORDER NUMBER
DS36 36 mm 	3.5	2	75 cm	blue		PP00DS36
	4	1	75 cm	blue		PP01DS36
DS40 40 mm 	3	2/0	75 cm	blue		PP20DS40
	3.5	0	75 cm	blue		PP00DS40
	4	1	90 cm	blue		PP01DS40

Reverse cutting needle – straight						
NEEDLE	EP	USP	LENGTH	SUTURE COLOUR	OTHER INFORMATION	ORDER NUMBER
GS60 60 mm 	2	3/0	75 cm	blue		PP30GS60
	3	2/0	75 cm	blue		PP20GS60

LUXAMID

NONABSORBABLE SURGICAL SUTURE

NYLON BLUE

TYPE	Monofilament
COMPOSITION	Extrusion of polyamide 6.0 or 6.6
COATING	None
TISSUE REACTION	Minimal
ABSORPTION	Non absorbable, gradually encapsulated by connective tissue. The thread mass diminishes, approximately 10% a year by rupture of chemical links (hydrolytic action).
SUTURE COLOUR	Blue, black
RANGE	USP 10/0 - 3&4
INDICATIONS	Luxamid sutures are intended for use in general soft tissue closing and/or ligation; especially in Luxamid blue: skin closure, general surgery, plastic surgery, gastrointestinal surgery, gynaecology, obstetrics, orthopaedics Luxamid black: plastic surgery, ophthalmic surgery





Round bodied needle – taper point – 1/2 circle

NEEDLE			EP	USP	LENGTH	SUTURE COLOUR	OTHER INFORMATION	ORDER NUMBER
HR16	16 mm		1.5	4/0	75 cm	blue		NY40HR16
HR20	20 mm		2	3/0	75 cm	blue		NY30HR20
HR26	26 mm		3	2/0	75 cm	blue		NY20HR26
			3.5	0	75 cm	blue		NY00HR26
HR30	30 mm		2	3/0	75 cm	blue		NY30HR30
			3	2/0	75 cm	blue		NY20HR30
HR40	40 mm		3.5	0	75 cm	blue		NY00HR40
			4	1	75 cm	blue		NY01HR40
			5	2	75 cm	blue		NY02HR40



Round bodied needle – taper point – 1/2 circle – loop

NEEDLE			EP	USP	LENGTH	SUTURE COLOUR	OTHER INFORMATION	ORDER NUMBER
HR40	40 mm		3.5	0	150 cm	blue	loop	NY00HR40-L
			4	1	150 cm	blue	loop	NY01HR40-L
HR50	50 mm		4	1	150 cm	blue	loop	NY01HR50-L



Reverse cutting needle – 3/8 circle

NEEDLE			EP	USP	LENGTH	SUTURE COLOUR	OTHER INFORMATION	ORDER NUMBER
DS12	12 mm		0.5	7/0	75 cm	black		NY70DS12
			0.7	6/0	75 cm	black		NY60DS12
			1	5/0	75 cm	black		NY50DS12
			1.5	4/0	75 cm	blue		NY40DS12
DS16	16 mm		0.7	6/0	45 cm	black		NY60DS16-45
			0.7	6/0	75 cm	black		NY60DS16
			1	5/0	45 cm	black		NY50DS16-45
			1	5/0	75 cm	black		NY50DS16
			1.5	4/0	45 cm	blue		NY40DS16-45
			1.5	4/0	75 cm	blue		NY40DS16
			2	3/0	75 cm	blue		NY30DS16

Reverse cutting needle – 3/8 circle						
NEEDLE	EP	USP	LENGTH	SUTURE COLOUR	OTHER INFORMATION	ORDER NUMBER
DS19 19 mm 	1	5/0	75 cm	black		NY50DS19
	1.5	4/0	45 cm	blue		NY40DS19-45
	1.5	4/0	75 cm	blue		NY40DS19
	2	3/0	75 cm	blue		NY30DS19
	3	2/0	75 cm	blue		NY20DS19
DS24 24 mm 	1.5	4/0	45 cm	blue		NY40DS24-45
	1.5	4/0	75 cm	blue		NY40DS24
	2	3/0	45 cm	blue		NY30DS24-45
	2	3/0	75 cm	blue		NY30DS24
	3	2/0	75 cm	blue		NY20DS24
	3.5	0	75 cm	blue		NY00DS24
DS26 26 mm 	2	3/0	75 cm	blue		NY30DS26
	3	2/0	75 cm	blue		NY20DS26
DS30 30 mm 	2	3/0	75 cm	blue		NY30DS30
	3	2/0	75 cm	blue		NY20DS30
	3.5	0	75 cm	blue		NY00DS30
	4	1	75 cm	blue		NY01DS30
DS36 36 mm 	2	3/0	75 cm	blue		NY30DS36
	3	2/0	75 cm	blue		NY20DS36
	3.5	0	75 cm	blue		NY00DS36
	4	1	75 cm	blue		NY01DS36
DS40 40 mm 	3	2/0	75 cm	blue		NY20DS40
	3	2/0	90 cm	blue		NY20DS40-90
	3.5	0	75 cm	blue		NY00DS40
	4	1	75 cm	blue		NY01DS40
	5	2	90 cm	blue		NY02DS40
DS90 90 mm 	3	2/0	100 cm	blue		NY20DS90
	4	1	90 cm	blue		NY01DS90

Reverse cutting needle – 1/2 circle						
NEEDLE	EP	USP	LENGTH	SUTURE COLOUR	OTHER INFORMATION	ORDER NUMBER
HS16 16 mm 	1	5/0	75 cm	black		NY50HS16
HS17 17 mm 	1.5	4/0	75 cm	blue		NY40HS17

Sutures in cassettes



EP	USP	LENGTH	SUTURE COLOUR	OTHER INFORMATION	ORDER NUMBER
1.5	4/0	100	blue		NY40100
2	3/0	100	blue		NY30100
3	2/0	100	blue		NY20100
3.5	0	100	blue		NY00100
4	1	100	blue		NY01100
5	2	75	blue		NY0275
6	3&4	50	blue		NY0350

PTFE

NONABSORBABLE SURGICAL SUTURE

POLYTETRAFLUOROETHYLENE MONOFILAMENT

TYPE	Monofilament
COMPOSITION	Polytetrafluoroethylene (PTFE)
COATING	None
TISSUE REACTION	Minimal
ABSORPTION	Non absorbable
SUTURE COLOUR	White
RANGE	USP 6/0 - 3/0
INDICATIONS	PTFE sutures are intended for use in general soft tissue closing and/or ligation; especially in skin closure, general surgery, plastic surgery, gastrointestinal surgery, gynaecology, obstetrics, orthopaedics, ophthalmic surgery and oral surgery.





Round bodied needle – taper point – 1/2 circle

NEEDLE			EP	USP	LENGTH	SUTURE COLOUR	OTHER INFORMATION	ORDER NUMBER
HR16	16 mm		0.7	6/0	60 cm	white	without CE	PT60HR16
			1	5/0	60 cm	white	without CE	PT50HR16
			1.5	4/0	60 cm	white	without CE	PT40HR16
HR20	20 mm		1	5/0	60 cm	white	without CE	PT50HR20
			1.5	4/0	60 cm	white	without CE	PT40HR20
			2	3/0	60 cm	white	without CE	PT30HR20



Reverse cutting needle – 3/8 circle

NEEDLE			EP	USP	LENGTH	SUTURE COLOUR	OTHER INFORMATION	ORDER NUMBER
DS12	12 mm		0.7	6/0	60 cm	white	without CE	PT60DS12
DS16	16 mm		0.7	6/0	60 cm	white	without CE	PT60DS16
			1	5/0	60 cm	white	without CE	PT50DS16
			1.5	4/0	60 cm	white	without CE	PT40DS16
DS19	19 mm		1	5/0	60 cm	white	without CE	PT50DS19
			1.5	4/0	60 cm	white	without CE	PT40DS19
			2	3/0	60 cm	white	without CE	PT30DS19

NONABSORBABLE SURGICAL SUTURE

SUPRAMID

POLYAMID BLACK

TYPE	Twisted and coated multifilament, sizes 5/0 and 6/0 monofilament
COMPOSITION	Polyamide 6.6 strands enclosed in a polyamide 6 sheet
COATING	None
TISSUE REACTION	Minimal
ABSORPTION	Non absorbable, gradually encapsulated by connective tissue
SUTURE COLOUR	Black
RANGE	USP 6/0 - 7
INDICATIONS	Supramid sutures are intended for use in general soft tissue closing and/or ligation; especially in skin closure





Round bodied needle – taper point – 1/2 circle

NEEDLE		EP	USP	LENGTH	SUTURE COLOUR	OTHER INFORMATION	ORDER NUMBER
HR16	16 mm		1	5/0	75 cm	black	SB50HR16
			1.5	4/0	75 cm	black	SB40HR16
HR20	20 mm		1.5	4/0	75 cm	black	SB40HR20



Round bodied needle – taper point – straight

NEEDLE		EP	USP	LENGTH	SUTURE COLOUR	OTHER INFORMATION	ORDER NUMBER
GR22	22 mm		1.5	4/0	50 cm	black	SB40GR22



Round bodied needle – tapercut point – 3/8 circle

NEEDLE		EP	USP	LENGTH	SUTURE COLOUR	OTHER INFORMATION	ORDER NUMBER
DRT18	18 mm		1.5	4/0	75 cm	black	SB40DRT18
			2	3/0	75 cm	black	SB30DRT18



Round bodied needle – tapercut point – 1/2 circle

NEEDLE		EP	USP	LENGTH	SUTURE COLOUR	OTHER INFORMATION	ORDER NUMBER
HRT16	16 mm		1	5/0	75 cm	black	SB50HRT16
HRT18	18 mm		1.5	4/0	75 cm	black	SB40HRT18
			2	3/0	75 cm	black	SB30HRT18

 Reverse cutting needle – 3/8 circle								
NEEDLE			EP	USP	LENGTH	SUTURE COLOUR	OTHER INFORMATION	ORDER NUMBER
DS12	12 mm		0.7	6/0	75 cm	black		SB60DS12
DS16	16 mm		0.7	6/0	75 cm	black		SB60DS16
			1	5/0	45 cm	black		SB50DS16-45
			1	5/0	75 cm	black		SB50DS16
			1.5	4/0	75 cm	black		SB40DS16-75
DS16	16 mm PREMIUM NEEDLE		1	5/0	75 cm	black		SB50DS16-P
DS19	19 mm		1	5/0	75 cm	black		SB50DS19
			1.5	4/0	45 cm	black		SB40DS19-45
			1.5	4/0	75 cm	black		SB40DS19
			2	3/0	75 cm	black		SB30DS19
DS24	24 mm		3	2/0	75 cm	black		SB20DS24

 Reverse cutting needle – 1/2 circle								
NEEDLE			EP	USP	LENGTH	SUTURE COLOUR	OTHER INFORMATION	ORDER NUMBER
HS16	16 mm		1	5/0	75 cm	black		SB50HS16
HS17	17 mm		1.5	4/0	75 cm	black		SB40HS17
HS20	20 mm		1	5/0	75 cm	black		SB50HS20
			1.5	4/0	75 cm	black		SB40HS20
			2	3/0	75 cm	black		SB30HS20
HS22	22 mm		1.5	4/0	75 cm	black		SB40HS22

Sutures in cassettes



EP	USP	LENGTH	SUTURE COLOUR	OTHER INFORMATION	ORDER NUMBER
1.5	4/0	100	black		SB40100
2	3/0	100	black		SB30100
3	2/0	100	black		SB20100
3.5	0	100	black		SB00100
4	1	100	black		SB01100
5	2	75	black		SB0275
6	3&4	50	black		SB0350
7	5	40	black		SB0540
8	6	40	black		SB0640
9	7	30	black		SB0730

NONABSORBABLE SURGICAL SUTURE

SUPRAMID

POLYAMID WHITE

TYPE	Twisted and coated multifilament
COMPOSITION	Polyamide 6.6 strands enclosed in a polyamide 6 sheet
COATING	None
TISSUE REACTION	Minimal
ABSORPTION	Non absorbable, gradually encapsulated by connective tissue
SUTURE COLOUR	White
RANGE	USP 4/0 - 8
INDICATIONS	Supramid sutures are intended for use in general soft tissue closing and/or ligation; especially in skin closure

Sutures in cassettes



EP	USP	LENGTH	SUTURE COLOUR	OTHER INFORMATION	ORDER NUMBER
1.5	4/0	100	white		SW40100
2	3/0	100	white		SW30100
3	2/0	100	white		SW20100
3.5	0	100	white		SW00100
4	1	100	white		SW01100
5	2	75	white		SW0275
6	3&4	50	white		SW0350
7	5	40	white		SW0540
8	6	40	white		SW0640
9	7	30	white		SW0730
10	8	30	white		SW0830

in solution

1.5	4/0	100	white		SWL40100
2	3/0	100	white		SWL30100
2.5	2/0	100	white		SWL20100
3	2/0	100	white		SWL20100
3.5	0	100	white		SWL00100
4	1	100	white		SWL01100
5	2	75	white		SWL0275
6	3&4	50	white		SWL0350
7	5	40	white		SWL0540
8	6	40	white		SWL0640
9	7	30	white		SWL0730
10	8	30	white		SWL0830

NONABSORBABLE SURGICAL SUTURE

LUXPET

POLYESTER BRAIDED

TYPE	Braided multifilament
COMPOSITION	Polyester - a polymer of polyethylene terephthalate
COATING	Silicone
TISSUE REACTION	Minimal
ABSORPTION	Non absorbable, gradually encapsulated by connective tissue
SUTURE COLOUR	Green White
RANGE	USP 5/0 - 8
INDICATIONS	Luxpet sutures are intended for use in general soft tissue closing and/ or ligation; especially in ophthalmic surgery, general surgery, skin closure, orthopaedics, gastrointestinal surgery, obstetrics



Without needles

NEEDLE	EP	USP	LENGTH	SUTURE COLOUR	OTHER INFORMATION	ORDER NUMBER
	3.5	0	150 cm	green		PE00WN150
	4	1	150 cm	green		PE01WN150
	5	2	150 cm	green		PE02WN150

Without needles – ligapack

NEEDLE	EP	USP	LENGTH	SUTURE COLOUR	OTHER INFORMATION	ORDER NUMBER
	2	3/0	250 cm	green		PE30WN250
	3	2/0	250 cm	green		PE20WN250



Round bodied needle – taper point – 1/2 circle

NEEDLE	EP	USP	LENGTH	SUTURE COLOUR	OTHER INFORMATION	ORDER NUMBER
HR16 16 mm 	1.5	4/0	75 cm	green		PE40HR16
HR20 20 mm 	2	3/0	75 cm	green		PE30HR20
HR26 26 mm 	2	3/0	75 cm	green		PE30HR26
	3	2/0	75 cm	green		PE20HR26
HR30 30 mm 	3	2/0	75 cm	green		PE20HR30
	3.5	0	75 cm	green		PE00HR30
HR36 36 mm 	3.5	0	75 cm	green		PE00HR36
	4	1	75 cm	green		PE01HR36
HR40 40 mm 	4	1	75 cm	green		PE01HR40
	5	2	75 cm	green		PE02HR40



Round bodied needle – taper point – 1/2 circle – double armed needle

NEEDLE	EP	USP	LENGTH	SUTURE COLOUR	OTHER INFORMATION	ORDER NUMBER
HR25 25 mm 	3	2/0	75 cm	green		PE20HR25-D



Round bodied needle – tapercut point – 1/2 circle

NEEDLE	EP	USP	LENGTH	SUTURE COLOUR	OTHER INFORMATION	ORDER NUMBER
HRT18 18 mm 	1.5	4/0	75 cm	green		PE40HRT18
	2	3/0	75 cm	green		PE30HRT18



Round bodied needle – tapercut point – 1/2 circle – double armed needle

NEEDLE	EP	USP	LENGTH	SUTURE COLOUR	OTHER INFORMATION	ORDER NUMBER
HRT17 17 mm 	3	2/0	90 cm	green		PE20HRT17-D



Round bodied needle – tapercut point – 1/2 circle – double armed needle

NEEDLE	EP	USP	LENGTH	SUTURE COLOUR	OTHER INFORMATION	ORDER NUMBER
HRT17 17 mm 	3	2/0	90 cm	white green	8 sutures / foil 5 foils / box with pledget 6 x 3 mm oval	PE20HRT17-DP
HRT25 25 mm 	3	2/0	90 cm	white green	8 sutures / foil 5 foils / box with pledget 6 x 3 mm oval	PE20HRT25-DP



Reverse cutting needle – 3/8 circle

NEEDLE	EP	USP	LENGTH	SUTURE COLOUR	OTHER INFORMATION	ORDER NUMBER
DS16 16 mm 	1	5/0	75 cm	green		PE50DS16
	1.5	4/0	75 cm	green		PE40DS16
DS19 19 mm 	1.5	4/0	75 cm	green		PE40DS19
	2	3/0	75 cm	green		PE30DS19
DS24 24 mm 	2	3/0	75 cm	green		PE30DS24
	3	2/0	75 cm	green		PE20DS24
DS30 30 mm 	3.5	0	75 cm	green		PE00DS30

Sutures in cassettes



EP	USP	LENGTH	SUTURE COLOUR	OTHER INFORMATION	ORDER NUMBER
1.5	4/0	100 m	green		PEG40100
2	3/0	100 m	green		PEG30100
3	2/0	100 m	green		PEG20100
3.5	0	100 m	green		PEG00100
4	1	100 m	green		PEG01100
5	2	75 m	green		PEG0275
6	3&4	50 m	green		PEG0350
7	5	40 m	green		PEG0540
8	6	40 m	green		PEG0640
3	2/0	100 m	white		PEW20100
3.5	0	100 m	white		PEW00100
4	1	100 m	white		PEW01100
5	2	100 m	white		PEW02100
6	3&4	80 m	white		PEW0380
7	5	60 m	white		PEW0560
8	6	50 m	white		PEW0650
9	7	40 m	white		PEW0740
10	8	30 m	white		PEW0830

NONABSORBABLE SURGICAL SUTURE

SILK BRAIDED

BLACK

TYPE	Braided multifilament
COMPOSITION	Braided fibres from the cocoon of the silkworm
COATING	Silicone
TISSUE REACTION	Moderate
ABSORPTION	Silk suture elicits an initial inflammatory reaction in tissues, which is followed by gradual encapsulation of the suture by fibrous connective tissues
SUTURE COLOUR	Black
RANGE	USP 8/0 - 8
INDICATIONS	Silk sutures are intended for use in general soft tissue closing and/or ligation; especially in general surgery, skin closure, ophthalmic surgery, oral surgery, gastrointestinal surgery, gynaecology, obstetrics



Without needles

NEEDLE	EP	USP	LENGTH	SUTURE COLOUR	OTHER INFORMATION	ORDER NUMBER
	2	3/0	150 cm	black		SI30WN150
	3	2/0	150 cm	black		SI20WN150
	3.5	0	150 cm	black		SI00WN150
	4	1	150 cm	black		SI01WN150
	2	3/0	13 x 60 cm	black		SI30WN60-13
	3	2/0	13 x 60 cm	black		SI20WN60-13
	3.5	0	13 x 60 cm	black		SI00WN60-13
	4	1	13 x 60 cm	black		SI01WN60-13
	3	2/0	12 x 75 cm	black		SI20WN75-12
	3.5	0	12 x 75 cm	black		SI00WN75-12
	4	1	12 x 75 cm	black		SI01WN75-12
	5	2	12 x 75 cm	black		SI02WN75-12

Without needles – ligapack

NEEDLE	EP	USP	LENGTH	SUTURE COLOUR	OTHER INFORMATION	ORDER NUMBER
	3	2/0	250 cm	black		SI20WN250
	4	1	250 cm	black		SI01WN250
	5	0	250 cm	black		SI02WN250



Round bodied needle – taper point – 1/2 circle

NEEDLE	EP	USP	LENGTH	SUTURE COLOUR	OTHER INFORMATION	ORDER NUMBER
HR13 13 mm 	1	5/0	45 cm	black		SI50HR13
HR16 16 mm 	1	5/0	75 cm	black		SI50HR16
	1.5	4/0	75 cm	black		SI40HR16
HR20 20 mm 	1	5/0	75 cm	black		SI50HR20
	1.5	4/0	75 cm	black		SI40HR20
	2	3/0	75 cm	black		SI30HR20
	3	2/0	75 cm	black		SI20HR20
HR26 26 mm 	1.5	4/0	75 cm	black		SI40HR26
	2	3/0	75 cm	black		SI30HR26
	3	2/0	75 cm	black		SI20HR26
	3.5	0	75 cm	black		SI00HR26
	4	1	75 cm	black		SI01HR26



Round bodied needle – taper point – 1/2 circle

NEEDLE		EP	USP	LENGTH	SUTURE COLOUR	OTHER INFORMATION	ORDER NUMBER
HR30 30 mm		2	3/0	75 cm	black		SI30HR30
		3	2/0	75 cm	black		SI20HR30
		3.5	0	75 cm	black		SI00HR30
		4	1	75 cm	black		SI01HR30
HR37 37 mm		3.5	0	75 cm	black		SI00HR37
		4	1	75 cm	black		SI01HR37
		5	2	75 cm	black		SI02HR37
HR40 40 mm		3.5	0	75 cm	black		SI00HR40
		4	1	75 cm	black		SI01HR40
		5	2	75 cm	black		SI02HR40



Round bodied needle – taper point – straight

NEEDLE		EP	USP	LENGTH	SUTURE COLOUR	OTHER INFORMATION	ORDER NUMBER
GR22 22 mm		1.5	4/0	50 cm	black		SI40GR22



Round bodied needle – tapercut point – 3/8 circle

NEEDLE		EP	USP	LENGTH	SUTURE COLOUR	OTHER INFORMATION	ORDER NUMBER
DRT18 18 mm		1.5	4/0	75 cm	black		SI40DRT18
		2	3/0	75 cm	black		SI30DRT18



Round bodied needle – tapercut point – 1/2 circle

NEEDLE		EP	USP	LENGTH	SUTURE COLOUR	OTHER INFORMATION	ORDER NUMBER
HRT18 18 mm		1.5	4/0	75 cm	black		SI40HRT18
		2	3/0	75 cm	black		SI30HRT18



Reverse cutting needle – 3/8 circle

NEEDLE			EP	USP	LENGTH	SUTURE COLOUR	OTHER INFORMATION	ORDER NUMBER
DS12	12 mm		0.7	6/0	45 cm	black		SI60DS12
			1	5/0	45 cm	black		SI50DS12
DS16	16 mm		0.7	6/0	75 cm	black		SI60DS16
			1	5/0	75 cm	black		SI50DS16
			1.5	4/0	75 cm	black		SI40DS16
			2	3/0	75 cm	black		SI30DS16
DS19	19 mm		1	5/0	75 cm	black		SI50DS19
			1.5	4/0	45 cm	black		SI40DS19-45
			1.5	4/0	75 cm	black		SI40DS19
			2	3/0	75 cm	black		SI30DS19
			3	2/0	75 cm	black		SI20DS19
DS24	24 mm		2	3/0	75 cm	black		SI30DS24
			3	2/0	75 cm	black		SI20DS24
			3.5	0	75 cm	black		SI00DS24
DS30	30 mm		1.5	4/0	75 cm	black		SI40DS30
			2	3/0	75 cm	black		SI30DS30
			3	2/0	75 cm	black		SI20DS30
			3.5	0	75 cm	black		SI00DS30
			4	1	75 cm	black		SI01DS30
DS37	37 mm		3	2/0	75 cm	black		SI20DS37
			3.5	0	75 cm	black		SI00DS37
			4	1	75 cm	black		SI01DS37
			5	2	75 cm	black		SI02DS37
DS45	45 mm		3.5	0	75 cm	black		SI00DS45
			4	1	75 cm	black		SI01DS45
			5	2	75 cm	black		SI02DS45



Reverse cutting needle – 1/2 circle

NEEDLE			EP	USP	LENGTH	SUTURE COLOUR	OTHER INFORMATION	ORDER NUMBER
HS17	17 mm		1.5	4/0	75 cm	black		SI40HS17
			2	3/0	75 cm	black		SI30HS17
HS20	20 mm		1.5	4/0	75 cm	black		SI40HS20
			2	3/0	75 cm	black		SI30HS20
HS22	22 mm		1.5	4/0	75 cm	black		SI40HS22
			2	3/0	75 cm	black		SI30HS22
			3	2/0	75 cm	black		SI20HS22
HS26	26 mm		2	3/0	75 cm	black		SI30HS26
			3	2/0	75 cm	black		SI20HS26

Sutures in cassettes

	EP	USP	LENGTH	SUTURE COLOUR	OTHER INFORMATION	ORDER NUMBER
	1	5/0	100 m	black		SI50100
	1.5	4/0	100 m	black		SI40100
	2	3/0	100 m	black		SI30100
	3	2/0	100 m	black		SI20100
	3.5	0	100 m	black		SI00100
	4	1	100 m	black		SI01100
	5	2	100 m	black		SI02100
	6	3&4	80 m	black		SI0380
	7	5	60 m	black		SI0560
	8	6	50 m	black		SI0650
	9	7	40 m	black		SI0740
	10	8	30 m	black		SI0830

STEEL

NONABSORBABLE SURGICAL SUTURE

STAINLESS STEEL MONOFILAMENT

TYPE	Monofilament
COMPOSITION	Iron, nickel and chromium alloy
COATING	None
TISSUE REACTION	Minimal
ABSORPTION	Non absorbable
SUTURE COLOUR	Metal
RANGE	USP 5/0 - 7
INDICATIONS	Steel sutures are intended for use in closure of sternum, orthopaedics, abdominal wall closure and hernia repair





Round bodied needle – tapercut point – 1/2 circle								
NEEDLE			EP	USP	LENGTH	SUTURE COLOUR	OTHER INFORMATION	ORDER NUMBER
HRT40	40 mm		6	3 & 4	4 x 45 cm	metal		ST03HRT40
HRT48	48 mm		7	5	4 x 45 cm	metal		ST05HRT48
			9	7	4 x 45 cm	metal		ST07HRT48



Reverse cutting needle – 1/2 circle								
NEEDLE		EP	USP	LENGTH	SUTURE COLOUR	OTHER INFORMATION	ORDER NUMBER	
HS40	40 mm		4	1	4 x 45 cm	metal		ST01HS40
			5	2	4 x 45 cm	metal		ST02HS40
HS48	48 mm		6	3 & 4	4 x 45 cm	metal		ST03HS48
HS50	50 mm		7	5	4 x 45 cm	metal		ST05HS50
			8	6	4 x 45 cm	metal		ST06HS50
HS58	58 mm		7	5	4 x 45 cm	metal		ST05HS58
			8	6	4 x 45 cm	metal		ST06HS58



Reverse cutting needle – 1/2 circle – PATELLA SET						
NEEDLE	EP	USP	LENGTH	SUTURE COLOUR	OTHER INFORMATION	ORDER NUMBER
NH120 120 mm	9	7	60 cm	metal	PATELLA SET	ST07HS120

Sutures in cassettes

	EP	USP	LENGTH	SUTURE COLOUR	OTHER INFORMATION	ORDER NUMBER
	1	5/0	50 m	metal		ST5050
	1.5	4/0	50 m	metal		ST4050
	2	3/0	50 m	metal		ST3050
	3	2/0	50 m	metal		ST2050
	3.5	0	50 m	metal		ST0050
	4	1	50 m	metal		ST0150

OPHTHALMIC + MICRO SUTURES



SYNTHETIC ABSORBABLE SURGICAL SUTURE
LUXCRYL 910
POLYGLACTINE 910 BRAIDED



Side cutting spatula – 1/2 circle – double needle

NEEDLE				EP	USP	LENGTH	SUTURE COLOUR	OTHER INFORMATION	ORDER NUMBER
HSP060	6 mm	200 µm		0.4	8/0	45 cm	violet		PL80HSP060-D
HSP080	8 mm	200 µm		0.5	7/0	45 cm	violet		PL70HSP080-D
HSP080	8 mm	300 µm		0.7	6/0	45 cm	violet		PL60HSP080-D



Side cutting spatula – 3/8 circle – double needle

NEEDLE				EP	USP	LENGTH	SUTURE COLOUR	OTHER INFORMATION	ORDER NUMBER
DSP060	6 mm	200 µm		0.4	8/0	45 cm	violet		PL80DSP060-D
				0.5	7/0	45 cm	violet		PL70DSP060-D
DSP065	6.5 mm	200 µm		0.4	8/0	45 cm	violet		PL80DSP065-D
				0.5	7/0	45 cm	violet		PL70DSP065-D
DSP065	6.5 mm	250 µm		0.7	6/0	45 cm	violet		PL60DSP065-D



Side cutting spatula – 1/4 circle – double needle

NEEDLE				EP	USP	LENGTH	SUTURE COLOUR	OTHER INFORMATION	ORDER NUMBER
VSP065	6.5 mm	300 µm		0.7	6/0	45 cm	violet		PL60VSP065-D
VSP080	8 mm	300 µm		0.7	6/0	45 cm	violet		PL60VSP080-D
				0.7	6/0	20 cm	undyed		PL60VSP080-DU
VSP080	8 mm	350 µm		1	5/0	45 cm	violet		PL50VSP080-D

OPHTHALMIC + MICRO SUTURES



SYNTHETIC ABSORBABLE SURGICAL SUTURE
LUXCRYL PGA
POLYGLYCOLIC ACID BRAIDED



Side cutting spatula – 3/8 circle – double needle

NEEDLE				EP	USP	LENGTH	SUTURE COLOUR	OTHER INFORMATION	ORDER NUMBER
DSP060	6 mm	200 µm		0.4	8/0	45 cm	violet		PG80DSP060-D
				0.5	7/0	45 cm	violet		PG70DSP060-D
DSP062	6.2 mm	150 µm		0.2	10/0	30 cm	violet		PG10DSP062-D
DSP064	6.4 mm	200 µm		0.3	9/0	30 cm	violet		PG90DSP064-D
DSP065	6.5 mm	250 µm		0.7	6/0	45 cm	violet		PG60DSP065-D

 Side cutting spatula – 1/4 circle – double needle									
NEEDLE				EP	USP	LENGTH	SUTURE COLOUR	OTHER INFORMATION	ORDER NUMBER
VSP080	8 mm	300 µm		0.7	6/0	45 cm	violet		PG60VSP080-D
VSP080	8 mm	350 µm		1	5/0	45 cm	violet		PG50VSP080-D

OPHTHALMIC + MICRO SUTURES



Side cutting spatula – 3/8 circle – double needle									
NEEDLE				EP	USP	LENGTH	SUTURE COLOUR	OTHER INFORMATION	ORDER NUMBER
DSP065	6.5 mm	200 µm		0.4	8/0	45 cm	black		SI80DSP065-D
DSP065	6.5 mm	250 µm		0.5	7/0	45 cm	black		SI70DSP065-D
				0.7	6/0	45 cm	black		SI60DSP065-D
DSP080	8 mm	350 µm		1	5/0	30 cm	black		SI50DSP080-D

Side cutting spatula – 1/4 circle – double needle									
NEEDLE				EP	USP	LENGTH	SUTURE COLOUR	OTHER INFORMATION	ORDER NUMBER
VSP080	8 mm	250 µm		0.7	6/0	45 cm	black		SI60VSP080-D



Side cutting spatula – 3/8 circle – double needle									
NEEDLE				EP	USP	LENGTH	SUTURE COLOUR	OTHER INFORMATION	ORDER NUMBER
DSP065	6.5 mm	200 µm		0.4	8/0	45 cm	blue		SV80DSP065-D

OPHTHALMIC + MICRO SUTURES



Side cutting spatula – 1/2 circle – double needle

NEEDLE	EP	USP	LENGTH	SUTURE COLOUR	OTHER INFORMATION	ORDER NUMBER
HSP062 6.2 mm 150 µm 	0.2	10/0	30 cm	black		NY10HSP062-D



Side cutting spatula – 3/8 circle – double needle

NEEDLE	EP	USP	LENGTH	SUTURE COLOUR	OTHER INFORMATION	ORDER NUMBER
DSP060 6 mm 150 µm 	0.2	10/0	30 cm	black		NY10DSP060-D
DSP062 6.2 mm 150 µm 	0.2	10/0	30 cm	black		NY10DSP062-D
	0.3	9/0	30 cm	black		NY90DSP062-D
	0.4	8/0	30 cm	black		NY80DSP062-D
DSP062 6.2 mm 200 µm 	0.5	7/0	30 cm	black		NY70DSP062-D
DSP064 6.4 mm 220 µm 	0.2	10/0	30 cm	black		NY10DSP064-D
	0.2	10/0	38 cm	black		NY10DSP064-D-38
DSP080 8 mm 300 µm 	0.7	6/0	45 cm	black		NY60DSP080-D



Side cutting spatula – 3/8 circle

NEEDLE	EP	USP	LENGTH	SUTURE COLOUR	OTHER INFORMATION	ORDER NUMBER
DSP062 6.2 mm 150 µm 	0.2	10/0	15 cm	black		NY10DSP062



Round bodied micro needle – taper point – 3/8 circle

NEEDLE	EP	USP	LENGTH	SUTURE COLOUR	OTHER INFORMATION	ORDER NUMBER
MDR041 4.1 mm 100 µm 	0,1	11/0	10 cm	black		NY11MDR041-100
MDR051 5.1 mm 100 µm 	0.2	10/0	15 cm	black		NY10MDR051-100
	0.3	9/0	15 cm	black		NY90MDR051-100
MDR051 5.1 mm 150 µm 	0.4	8/0	15 cm	black		NY80MDR051-150
MDR064 6.4 mm 150 µm 	0.4	8/0	15 cm	black		NY80MDR064-150

OPHTHALMIC + MICRO SUTURES



NONABSORBABLE SURGICAL SUTURE

POLYESTER BRAIDED



Side cutting spatula – 1/4 circle – double needle

NEEDLE				EP	USP	LENGTH	SUTURE COLOUR	OTHER INFORMATION	ORDER NUMBER
VSP080	8 mm	350 µm		1	5/0	45 cm	white		PE50VSP080-D



NONABSORBABLE SURGICAL SUTURE

POLYPROPYLENE MONOFILAMENT



Side cutting spatula – 3/8 circle – double needle

NEEDLE				EP	USP	LENGTH	SUTURE COLOUR	OTHER INFORMATION	ORDER NUMBER
DSP064	6.4 mm	200 µm		0.2	10/0	30 cm	blue		PP10DSP064-D
				0,3	9/0	30 cm	blue		PP90DSP064-D



Side cutting spatula – straight

NEEDLE				EP	USP	LENGTH	SUTURE COLOUR	OTHER INFORMATION	ORDER NUMBER
GSP16	16.15 mm	150 µm		0.2	10/0	20 cm	blue		PP10GSP16-D

SURGICAL NEEDLES




1/2 CIRCLE

12 in one box

Reverse cutting needle

Round bodied needle



	Ordinary Eye GN	Spring Eye GF	Ordinary Eye GRN	Spring Eye GRF
90 mm N°00		1HS90-S		
80 mm N°1	1HS80-O	1HS80-S		1HR80-S
70 mm N°2	1SH70-O	1HS70-S		1HR70-S
65 mm N°3	1HS65-O			
60 mm N°4		1HS60-S		1HR60-S
55 mm N°6	1HS55-O	1HS55-S	1HR55-O	1HR55-S
50 mm N°7		1HS50-S	1HR50-O	1HR50-S
45 mm N°8		1HS45-S		1HR45-S
40 mm N°9	1HS40-O	1HS40-S	1HR40-O	1HR40-S
35 mm N°10		1HS35-S		1HR35-S
30 mm N°12		1HS30-S	1HR30-O	1HR30-S
25 mm N°13	1HS25-O	1HS25-S	1HR25-O	1HR25-S
20 mm N°15		1HS20-S		1HR20-S
16 mm N°16		1HS16-S		1HR16-S

SKI NEEDLES half curved

12 in one box

Reverse cutting needle

Ordinary Eye

Spring Eye



	Ordinary Eye	Spring Eye
80 mm N°1	1SKI80-O	
70 mm N°2	1SKI70-O	
65 mm N°3	1SKI65-O	
60 mm N°4	1SKI60-O	
57 mm N°5	1SKI57-O	
55 mm N°6	1SKI55-O	
50 mm N°7	1SKI50-O	
45 mm N°8	1SKI45-O	
40 mm N°9		1SKI40-S
35 mm N°10		1SKI35-S
30 mm N°12		1SKI30-S
25 mm N°13		1SKI25-S
20 mm N°15		1SKI20-S
16 mm N°16		1SKI16-S
12 mm N°20		1SKI12-S

SERPENTINE NEEDLES

12 in one box

Reverse cutting needle

Ordinary Eye



	Ordinary Eye
135 mm	1SERP135-O
145 mm	1SERP145-O

3/8 CIRCLE

12 in one box



Reverse cutting needle

Round bodied
taper point

needle
tapercut point

Ordinary
Eye

Spring
Eye

Ordinary
Eye

Spring
Eye

Ordinary
Eye

130 mm N°6-0	1DS130-O		1DR130-O		
110 mm N°4-0			1DR110-O		
100 mm N°000	1DS100-O	1DS100-S		1DR100-S	
90 mm N°00	1DS90-O	1DS90-S	1DR90-O	1DR90-S	1DRT90-O
90 mm LIGHT N°00	1DS90-OL		1DR90-OL		
85 mm N°0			1DR85-O		
80 mm N°1	1DS80-O	1DS80-S	1DR80-O		1DRT80-O
70 mm N°2	1DS70-O	1DS70-S	1DR70-O	1DR70-S	
65 mm N°3			1DR65-O		
60 mm N°4	1DS60-O	1DS60-S		1DR60-S	
57 mm N°5	1DS57-O		1DR57-O		
55 mm N°6		1DS55-S			
50 mm N°7		1DS50-S	1DR50-O	1DR50-S	
45 mm N°8	1DS45-O	1DS45-S	1DR45-O	1DR45-S	
40 mm N°9	1DS40-O	1DS40-S	1DR40-O	1DR40-S	
35 mm N°10		1DS35-S		1DR35-S	
30 mm N°12	1DS30-O	1DS30-S	1DR30-O	1DR30-S	
25 mm N°13	1DS25-O	1DS25-S		1DR25-S	
20 mm N°15		1DS20-S		1DR20-S	
16 mm N°16	1DS16-O	1DS16-S		1DR16-S	
14 mm N°17		1DS14-S			
12 mm N°20		1DS12-S			

STRAIGHT NEEDLES

12 in one box

Reverse cutting needle


Ordinary Eye

Spring Eye

75 mm N° 1		1GS75-S
60 mm N° 4	1GS60-O	1GS60-S
45 mm N° 8	1GS45-O	
40 mm N° 9		1GS40-S
35 mm N° 10		1GS36-S

POST MORTEM NEEDLE 3/8 circle

12 in one box



Reverse cutting needle

Ordinary Eye

150 mm N°6"	1DS150-O
130 mm N°5"	1DS130-O
100 mm N°4"	1DS100-O

OTHER PRODUCTS



SURGICAL BLADES



Surgical blades with plastic handle



Individual sterile foil packs
without handle 100 / box
with plastic handle 10 / box



COMPOSITION:
Carbon steel

STITCH CUTTER



LUSC

COMPOSITION:
Carbon steel

Individual sterile foil packs 100 / box



BONEWAX 2.5 g
SURGICAL HAEMOSTATIC
12 x 2.5 g **LUBW**

COMPOSITION:
Refined white beeswax, paraffin and isopropyl palmitate

INDICATIONS:
Used in the control of bleeding from bone surfaces by acting as a mechanical barrier.

SURGICAL MESHES

Polypropylene hernia mesh


15 x 15 cm	LUMP1515
30 x 30 cm	LUMP3030
6 x 11 cm	LUMP0611
10 x 15 cm	LUMP1015
7,5 x 15 cm	LUMP7515

COMPOSITION:

Monofilament polypropylene knitted into an elastic, durable mesh

INDICATIONS:

Reinforcement of the abdominal wall
- hernia
- eventration
via celioscopy or laparotomy



PILOMAT

Tablet crusher

Z045



Tablet Splitter

Z044



SKIN STAPLER

Skin stapler 35 W

LUSS

10 pc.



SUTURE RACK



Rack to join 125 mm x 72 mm x 143 mm
Z070 1 pc.

Drawer for rack
Z071 1 pc.

FLEX-BANDAGES

Complete range of flexible, self adhesive and comfortable bandages.
Suitable for all types of animals.

- Adheres only to itself, not to skin or clothing
- No pins or clips are required to hold the bandage in place
- Strong and pleasant to apply
- Easy to cut, no scissors needed
- Flexible substrate - allows the bandage to conform to awkward body contours
- Air permeable - allows the skin to breathe



2,5 cm x 4,5 m	mix with different colours	LUBAM25	60 pc.
5 cm x 4,5 m	mix with different colours	LUBAM5	36 pc.
7,5 cm x 4,5 m	mix with different colours	LUBAM75	24 pc.
10 cm x 4,5 m	mix with different colours	LUBAM10	18 pc.



2,5 cm x 4,5 m	special mix - paws and smileys	LUBAM25-Z	60 pc.
5 cm x 4,5 m	special mix - paws and smileys	LUBAM5-Z	36 pc.
7,5 cm x 4,5 m	special mix - paws and smileys	LUBAM75-Z	24 pc.
10 cm x 4,5 m	special mix - paws and smileys	LUBAM10-Z	18 pc.

LUXBOND – SKIN ADHESIVE



Luxbond Tissue Adhesive provides quality wound management. It polymerizes in seconds after contact with tissue and binds wound edges together so natural healing can occur.

- Good visibility due to blue colour
- Acts within seconds
- For additional suture support
- Protection against infections
- Stops Bleeding

Luxbond 2,5 ml	LUGV	1 pc.
Luxbond applicators	LUGVA	6 pc.

MAGNETS

magnet a cage green	LU080	12 pc.
magnet super fort	LU081	12 pc.
magnet super fort a cage green	LU081	12 pc.

PROTECTION

Body protection

- Disposable
- Long funnel-shaped raglan sleeves with tight forearms cuffs to provide a good seal
- Individually packed and folded gowns
- Long and roomy for kneeling



ECONOMIC 25 μ

180 cm green
170 cm green

24 pc. **ZP001**
24 pc. **ZP002**

STANDARD 30 μ

190 cm (165 cm) blue
180 cm (156 cm) blue
170 cm (140 cm) blue
180 cm (156 cm) blue sterile
170 cm (140 cm) blue sterile

24 pc. **ZP004**
24 pc. **ZP005**
24 pc. **ZP006**
1 pc. **ZP007**
1 pc. **ZP008**



Overall

PP-spun non-woven with cap, disposable suit, green

S	5 pc.	ZP050
M	5 pc.	ZP051
L	5 pc.	ZP052
XL	5 pc.	ZP053
XXL	5 pc.	ZP054



Overboots

OVERBOOTS 60 μ , 35 CM LONG AND 50 CM HIGH
Disposable plastic overboots with elastic bands

50 pc. **ZP030**





22, Gruuss-Strooss, L-9991 Weiswampach, Luxembourg
info@luxsutures.com www.luxsutures.com